

HIGH FREQUENCY **INNOVATION**
























Ajwa Surgical



















2020
CATALOGUE

Info@ajwasurgical.com

www.ajwasurgical.com

TITANIUM STANDARD	TEMPORARY						PERMANENT					
	Item No.		Blade length mm	Maximal opening mm	Closing force N g		Item No.		Blade length mm	Maximal opening mm	Closing force N g	
	FT240T		7.0	6.2	1.08	110	FT740T		7.0	6.2	1.47	150
	FT250T		9.0	7.0	0.88	90	FT750T		9.0	7.0	1.77	180
	FT260T		11.0	7.8	0.88	90	FT760T		11.0	7.8	1.77	180
	FT280T		15.0	9.2	0.88	90	FT780D		15.0	9.2	1.96	200
	FT292T		17.5	10.6	0.88	90	FT792D		17.5	10.6	1.96	200
	FT290T		20.0	11.4	1.08	110	FT790D		20.0	11.4	1.96	200
	FT242T		6.5	6.0	1.08	110	FT742T		6.5	6.0	1.47	150
	FT252T		8.3	6.8	0.88	90	FT752T		8.3	6.8	1.77	180
	FT262T		10.2	7.5	0.88	90	FT762T		10.2	7.5	1.77	180
	FT282T		13.7	8.7	1.08	110	FT782T		13.7	8.7	1.96	200
	FT286T		15.3	9.8	1.08	110	FT786T		15.3	9.8	1.77	180
	FT244T		5.4	5.5	1.08	110	FT744T		5.4	5.5	1.47	150
	FT254T		6.4	6.1	1.08	110	FT754T		6.4	6.1	1.77	180



TEMPORARY						PERMANENT					
Item No.		Blade length mm	Maximal opening mm	Closing force N g		Item No.		Blade length mm	Maximal opening mm	Closing force N g	
FT264T		8.0	6.5	1.08	110	FT764T		8.0	6.5	1.77	180
FT284T		10.3	7.4	1.08	110	FT784T		10.3	7.4	1.96	200
FT323T		7.4	5.4	1.08	110	FT823T		7.4	5.4	1.77	180
FT324T		8.4	5.8	1.08	110	FT824T		8.4	5.8	1.77	180
FT325T		10.5	6.2	1.08	110	FT825T		10.5	6.2	1.77	180
FT270T		5.0	8.0	1.08	110	FT770T		5.0	8.0	1.96	200
FT271T		5.0	8.0	1.08	110	FT771T		5.0	8.0	1.96	200
FT272T		5.0	8.0	1.08	110	FT772T		5.0	8.0	1.96	200

■ MR Conditional, 3-Tesla or less*

■ Color coding with blue spring (Standard clip) and golden or silver blades (temporary or permanent clip respectively)

■ Sterile packaging

* Please refer to the MRI Safety Information section of the Instructions for Use for MR-related safety information.















TEMPORARY						PERMANENT					
Item No.		Blade length mm	Maximal opening mm	Closing force N g		Item No.		Blade length mm	Maximal opening mm	Closing force N g	
FT322T		10.0	5.6	1.28	130	FT822T		10.0	5.6	1.96	200
FT350T		6.7	5.4	1.28	130	FT850T		6.7	5.4	1.96	200
FT351T		9.0	5.6	1.28	130	FT851T		9.0	5.6	1.96	200
FT330T		7.0	7.2	0.88	90	FT830T		7.0	7.2	1.96	200
FT332T		9.3	7.2	0.88	90	FT832T		9.3	7.2	1.96	200
FT333T		11.3	8.2	0.88	90	FT833T		11.3	8.2	1.96	200
FT246T		7.0	5.7	0.88	90	FT746T		7.0	5.7	1.96	200
FT273T		7.8	5.3	0.88	90	FT773T		7.8	5.5	1.77	180
FT261T		11.4	8.0	0.88	90	FT761T		11.4	8.0	1.96	200
FT247T		8.6	7.0	0.88	90	FT747T		8.6	7.0	1.77	180

- MR Conditional, 3-Tesla or less*
- Color coding with blue spring (Standard clip) and golden or silver blades (temporary or permanent clip respectively)
- Sterile packaging
- * Please refer to the MRI Safety Information section of the Instructions for Use for MR-related safety information.

TEMPORARY

PERMANENT

Item No.		Diameter of fenestration 3.5 mm				Item No.		Diameter of fenestration 3.5 mm			
		Blade length mm L1/L2	Maximal opening mm	Closing force N g				Blade length mm L1/L2	Maximal opening mm	Closing force N g	
FT097T		3/8.1	7.5	1.08	110	FT597T		3/8.1	7.5	1.47	150
FT098T		4/9.1	7.4	1.08	110	FT598T		4/9.1	7.4	1.47	150
FT099T		5/10.1	7.4	1.08	110	FT599T		5/10.1	7.4	1.47	150
FT100T		6/11.1	7.4	1.08	110	FT600T		6/11.1	7.4	1.47	150
FT110T		9/14.1	8.4	1.08	110	FT610T		9/14.1	8.4	1.47	150
FT120T		12/17	9.6	1.08	110	FT620T		12/17	9.6	1.77	180
FT102T		5/7.8	5.8	1.08	110	FT602T		5/7.8	5.8	1.47	150
FT112T		7.5/9.4	6.5	1.08	110	FT612T		7.5/9.4	6.5	1.47	150
FT122T		10/11	7.2	1.08	110	FT622T		10/11	7.2	1.77	180
FT104T		5/5.1	4.9	1.08	110	FT604T		5/5.1	4.9	1.47	150
FT114T		7.5/5.1	4.9	1.08	110	FT614T		7.5/5.1	4.9	1.47	150

STANDARD

MINI



FT501T

FT511T

with latch

FT505T

FT515T

without latch



Total length: 225 mm, 8 7/8 inch



FT502T


FT512T

with latch

FT506T

FT516T

without latch



Total length: 245 mm, 9 5/8 inch



- The Double Bayonet Clip Applier has an additional bayonet shape at the jaw.



- This allows for a better view of the clip.

■ Color coded is indicated on the handle (● = Standard; ● = Mini)

STANDARD

MINI



FT501T

FT511T

with latch

FT505T

FT515T

without latch



Total length: 225 mm, 8 7/8 inch



FT502T


FT512T

with latch

FT506T

FT516T

without latch



Total length: 245 mm, 9 5/8 inch



- The Double Bayonet Clip Applier has an additional bayonet shape at the jaw.



- This allows for a better view of the clip.

■ Color coded is indicated on the handle (● = Standard; ● = Mini)



Ajwa Surgical

Mohallah Jeevan Pura Kacha Shahabpura, Sialkot (51310) Pakistan.

Phone: 0092 (315) 7651238

E-mail: info@ajwasurgical.com

URL: [http:// www.ajwasurgical.com](http://www.ajwasurgical.com)